


SampliQ QuEChERS 買一送一活動

活動時間 2009/4/15 至 2009/9/30

Buy any Agilent SampliQ QuEChERS
Extraction Kit and Receive an
Agilent SampliQ QuEChERS
Dispersive SPE Kit FREE



Agilent SampliQ QuEChERS Kits

Save time and money with a simplified approach to preparing samples for pesticide analysis.

In 2002, USEPA scientists developed a groundbreaking method for simplifying the very time-consuming task of sample extraction. It's called QuEChERS (quEChERS = quick, easy, cheap, efficient, rugged and safe). With QuEChERS, you can prepare your samples for analysis with fewer, smaller vials and less time.

See how by more information about the method... and the way that Agilent SampliQ makes it easy for you to quickly start using QuEChERS to improve your productivity.



Our mission is your success.

products | applications | software | services



Agilent Technologies

歡
迎
來
電
洽
詢



Thank you for trying Agilent SampliQ QuEChERS Kits

Agilent SampliQ
High performance solid phase extraction products

Now you can save time with a simplified approach to preparing samples for pesticide analysis.

The QuEChERS method can help you save significant time in food sample preparation for pesticide analysis. With SampliQ QuEChERS kits, you have no weighing or spiking, making the method even easier and more consistent.



Get started now

Place any of the boxes shown to complete the QuEChERS method on any food or vegetable sample. To get a free video showing you how to prepare your samples for multi-residue, multi-residue pesticide analysis.

STEP ONE: SampliQ QuEChERS Selection Kits

QuEChERS Kits have 50 mL centrifuge tubes, 50 mL vials, seals, centrifuge products, containing pre-measured, pre-mixed salts as listed.

Description	Part No.	Part No.
Standard QuEChERS Extraction Kit, 5000 Medical 500-01	5000-0100	5000-0100
Standard QuEChERS Extraction Tube, 50 Medical 500-02	5000-0200	5000-0200
Agilent QuEChERS Method (see Informed Extraction Tubes, 50 gram)	5000-0300	5000-0300
Agilent QuEChERS Method (see Informed Extraction Tubes, 10 gram)	5000-0400	5000-0400
QuEChERS Informed Tubes (see Informed)	5000-0500	5000-0500

STEP TWO: SampliQ QuEChERS Dispersive SPE Kits

QuEChERS Kits have 50 mL centrifuge tubes, 50 mL vials, seals, centrifuge products, containing pre-measured, pre-mixed salts as listed.

Sample Type	Tube Size (mL)	Part No. (50)	Part No. (10)
Standard Kit	2.0 mL (100 µL)	5000-0601	5000-0602
Standard Kit	15 mL (150 µL)	5000-0603	5000-0604
Standard Kit	15 mL (150 µL)	5000-0605	5000-0606
Standard Kit	15 mL (150 µL)	5000-0607	5000-0608
Standard Kit	15 mL (150 µL)	5000-0609	5000-0610
Standard Kit	15 mL (150 µL)	5000-0611	5000-0612
Standard Kit	15 mL (150 µL)	5000-0613	5000-0614
Standard Kit	15 mL (150 µL)	5000-0615	5000-0616
Standard Kit	15 mL (150 µL)	5000-0617	5000-0618
Standard Kit	15 mL (150 µL)	5000-0619	5000-0620

See our QuEChERS website, and view videos of the QuEChERS method.

Take your QuEChERS method to the next level. Find out about special offers, and order kits at www.agilent.com/chem/quchers

Go to the bottom of product pages on the Agilent Technical Support at 1-800-352-3700 or the U.S. and Canada Toll-Free number, and www.agilent.com/chem/quchers

Our mission is your success.

Agilent Technologies



巨研科技全省耗材服務諮詢專線：

北部：新竹以北(含) (02)29457272

轉 304 施小姐

中部：苗栗以南、嘉義以北 (07)5507380

轉 216 林小姐

南部：台南以南(含)、高屏地區(07)5507380

轉 209 陳小姐

The new color coded manual syringes 活動

活動優惠: Buy 4 syringes get 25% discount on the order

活動時間 2009/4/15 至 2009/9/30

Agilent is changing the way you see syringes.



**Introducing
NEW Agilent Manual Syringes**

New, color-coded manual syringes let you identify critical volumes at a glance and help boost productivity along every step of your sample processing workflow.

Whether performing manual dilution, extraction, enhancement or sample preparation, deciphering tiny part numbers on syringe barrels is a waste of time in today's busy labs. That's why Agilent is pleased to introduce a new line of syringes that are color-coded based on their maximum volume capacity.



Each of the manual syringes has a distinct color bonded onto every glass barrel so you have the every time. The syringes also come with a convenient stand for safe storage and are delivered in easy-to-use, recyclable packaging.

Color Speaks Volumes!

Our measure is your success.

 **Agilent Technologies**

You'll never look at another syringe the same way.



The difference should be easy to see:

- Available in volumes from 0.5 μ L to 50 mL
- Distinctive color permanently bonded onto every glass barrel for fast, easy identification
- Employ an intuitive, easy-to-read, vertical orientated scale
- NIST traceable
- Provide accuracy \pm 1% for Plunger-in-Barrel syringes and \pm 2% for Plunger-in-Needle syringes
- Use innovative packaging that doubles as a bench top storage stand
- Use environmentally friendly recyclable packaging
- Lot tracking and part number bonded directly on syringe

Barrel Color	Volume Size
Light Blue	0.5 μ L
Light Red	1.0 μ L
Light Green	2.0 μ L
Light Yellow	5.0 μ L
Light Purple	10 μ L / 1.0 mL
Light Orange	25 μ L / 2.5 mL
Light Blue-Grey	50 μ L / 5.0 mL
Light Yellow-Green	100 μ L / 10.0 mL
Light Green	250 μ L / 25.0 mL
Light Blue	500 μ L / 50.0 mL

Agilent syringes can be ordered in a variety of flexible and plunger configurations to meet all of your sample preparation needs. Visit www.agilent.com/chem/syringes to learn more.

Agilent manual syringes are available in convenient singles and economical 5-packs and 10-packs for high volume needs. **One look and you'll see why new Agilent manual syringes are a vivid, accurate and brilliant addition to your lab.**

Download a copy of Agilent's Syringe Specifications and Selection Guide today!

With a broad selection of syringes for manual and auto injectors, Agilent has what you need for accurate and effective sampling.

Free product specifications, selection guides and handy tips on how to determine the right syringe for your application. Our 24 page booklet covers the complete line of Agilent manual and autosampler syringes. To download our free Syringe Selection Guide, visit:

www.agilent.com/chem/syringes

 **Agilent Technologies**

巨研科技全省耗材服務諮詢專線:

北部:新竹以北(含) (02)29457272 轉 304 施小姐

中部:苗栗以南、嘉義以北 (07)5507380 轉 216 林小姐

南部:台南以南(含)、高屏地區(07)5507380 轉 209 陳小姐

『ZORBAX Solvent Saver Columns 優惠促銷活動』

ACN 是一種非常廣泛應用在液相層析儀的溶劑，因為石化製造過程需要使用 ACN 合成纖維，現今 ACN 的供應全球化短缺，價格飆昇，為了減少 ACN 的使用量，在既有的 HPLC 儀器來看，主要朝管柱內徑改小(流量變小)、粒子改小(時間縮短)、長度改短的方向改善。

Agilent 為了因應減少 ACN 使用量，建議客戶使用 Solvent Saver 管柱/Rapid Resolution 管柱，幫助客戶減少 50% 的溶劑或更多的溶劑。

例如從 4.6x150mm, 5 μ m \rightarrow 3.0x150mm, 3.5 μ m 的來看:

- (1) 減少溶劑至少 50%，節省成本。
- (2) 更換孔徑 3.0mm 管柱不需更改系統(flow cell、tubing)。
- (3) 解析度提升、時間縮短，提升工作效率。

自即日起至 2009/10/15，凡購買 Solvent Saver (3.0mm ID) 或 Narrow

Bore (2.1mm ID) 的 ZORBAX Columns 可享有 30% 折扣，歡迎來電洽詢。



巨研科技全省耗材服務諮詢專線:

北部: 新竹以北(含) (02)29457272 轉 304 施小姐
中部: 苗栗以南、嘉義以北 (07)5507380 轉 216 林小姐
南部: 台南以南(含)、高屏地區(07)5507380 轉 209 陳小姐